



Sienna COVID-19 IgG/IgM Rapid Test

Test Details

Sienna COVID-19 IgG/IgM Rapid Test Cassette is designed to detect the presence of IgG and IgM antibodies to SARS-CoV-2 in blood, serum or plasma samples. This COVID-19 IgG/IgM Rapid Test is intended to aid in identifying individuals with an adaptive immune response to SARS-CoV-2, indicating a prior infection. The COVID-19 IgG/IgM Rapid Test Cassette is not to be used to diagnose acute SARS-CoV-2 infection.

Accurate Results

- Sensitivity: 93.3%
- Specificity: 98.8%
- Accuracy: 99.6%

Unmatched Convenience

- Visual Results: Easy to interpret.
- Fast Results: 10- 20 minute test time.

